


ALSEE

Glistening-Free Hydrophobic IOL



MONOFOCAL	Hydrophobic, Aspheric, Modified C
Optic	6.00 mm, Biconvex, Aspheric, Aberration Correction (-0,165 µm)
Overall Ø	13.00 mm, Modified C-Loop
Haptic Angle	0°
Recommended Ac. A Constant	118.0
Recommended Op. A Constant	SRK-T: 118.4 - SRK-II: 118.6
Diopter Power Range	From 0.0 D to +32.00 D (0.50 D increments)
Special Production	From -5.00 D to -0.50 D and +32.50 D to +45.00 D
Refractive Index	20°C /35°C 1.531 / 1.533 ± 0.002 @546 nm
Injection System	AlsaJET (1.8 mm and 2.2 mm)



ALSANZA
Medizintechnik und Pharma GmbH


www.alsanza.com

Hermann-Burkhardt-Straße 3 72793 Pfullingen Germany
info@alsanza.com



MONOFOCAL

ALSEE



**a glistening free hydrophobic
acrylic IOL packed in dry conditions**


Safe platform with a performing optic design

360° square edge and reliable design

Foldability - resistance to distortion

Premium aspheric optic concept

MORE INFO



ALSANZA
Medizintechnik und Pharma GmbH